

# helica Ochratoxin A ELISA

For Coffee, Spices, And More



[Catalog No. 961OCH01COF]

## Product Highlights

- Assay range 1-20ppb
- Limit of detection < 1ppb
- Easy extraction
- Cost-effective and user-friendly tests
- Results in 70 minutes
- Break-away wells allow low-volume testing
- Standards and conjugate preserved with proprietary stabilizers
- One year stability and reproducibility of results
- Easy to follow instructions
- Color-coded reagents in convenient, easy to use vials
- No confirmation by HPLC Required

## Mycotoxin ELISA KITS

| Catalog No.    | Product                                                                                               |
|----------------|-------------------------------------------------------------------------------------------------------|
| 961AFLM01M-96  | <b>Aflatoxin M1</b><br>for milk and dairy products                                                    |
| 941BAFL01B-96  | <b>Aflatoxin B1</b><br>for grains, cereal, nuts, & animal feed                                        |
| 941AFL01M-96   | <b>Total Aflatoxins</b><br>for grains, cereal, nuts, & animal feed                                    |
| 981AFL01LM-96  | <b>Total Aflatoxins Low Matrix</b><br>for grains, cereal, silage, nuts, spices,<br>& animal feed      |
| 991AFLM01U-96  | <b>Aflatoxin M1</b><br>For Urine                                                                      |
| 941DON01M-96   | <b>Deoxynivalenol (DON)</b><br>for grains, cereal, & animal feed                                      |
| 951FUM01C-96   | <b>Fumonisin</b><br>for maize (corn)                                                                  |
| 951ZEA01G-96   | <b>NEW Zearalenone</b><br>for corn, barley, oats, wheat,<br>rice, sorghum                             |
| 941OCH01M-96   | <b>Ochratoxin A</b><br>For grains, cereals,<br>animal feed- <b>EU Standards</b>                       |
| 981OCH01W-96   | <b>Ochratoxin A</b><br>Qualitative for wine, grape juice<br>& grape must                              |
| 981OCH01ALC-96 | <b>Ochratoxin A</b><br>Quantitative for alcoholic beverages,<br>wine, grape juice, grape must, & beer |
| 991OCH01MS-96  | <b>Ochratoxin A</b><br>for human and animal milk, serum,<br>& plasma                                  |
| 961OCH01COF-96 | <b>Ochratoxin A</b><br>for coffee, cocoa butter and spices                                            |



Helica Biosystems, Inc.  
1527 W. Alton Ave  
Santa Ana, CA 92704, USA

P: 714-578-7830

F: 714-578-7831

www.helica.com

info@helica.com